# DIAQUICK COVID-19 IgG/IgM Cassette

Whole Blood, Serum and Plasma

## REF

### Content

- Z20101CE
- 25 test cassettes, individually packed in foil pouches with a dropper and a desiccant (25x REF Z20101B)
  - 1 Buffer
  - 1 package insert

**For professional in vitro diagnostic use only.**

### GENERAL INFORMATION

- Method Lateral flow immunochromatographic assay  
 Storage 2-30 °C

### INTENDED USE

COVID-19 IgG/IgM Rapid Test Cassette is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies to 2019 Novel Coronavirus in human whole blood, serum or plasma. This test provides only preliminary test result. Therefore, any reactive specimen with the COVID-19 IgG/IgM Rapid Test Cassette must be confirmed with alternative testing method(s) and clinical findings.

### DIAGNOSTIC SIGNIFICANCE

Coronaviruses are enveloped RNA viruses that are common among humans, other mammals and birds and cause respiratory, enteric, hepatic and neurologic diseases. Seven coronavirus species are known to cause disease in humans. Four viruses- 229E, OC43, NL63 and HKU1 – are prevalent and typically cause common cold symptoms in immunocompetent individuals. The other three strains- severe acute respiratory syndrome coronavirus (SARS-CoV), Middle East respiratory syndrome Coronavirus (MERS-CoV) and 2019 Novel Coronavirus (COVID-19)- are of zoonotic origin and have been linked to sometimes fatal diseases. IgG and IgM antibodies against the 2019 Novel Coronavirus can be detected 2-3 weeks after exposure. IgG remains positive, but antibody levels decrease over time.

### TEST PRINCIPLE

The COVID-19 IgG/IgM Rapid Test Cassette is a lateral flow immunochromatographic assay. The test uses anti-human IgM antibody (test line IgM), anti-human IgG (test line IgG) and goat anti-rabbit IgG (control line C) immobilised on nitrocellulose strip. The burgundy colored conjugate pad contains colloidal gold conjugated to recombinant COVID-19 antigens conjugated with colloid gold (COVID-19 conjugates) and rabbit IgG-gold conjugates. When a sample followed by assay buffer is added to the sample well, IgM and/or IgG antibodies, if present, bind to COVID-19 conjugates, thereby forming antigen-antibody complexes. This complex migrates through the capillary action across the nitrocellulose membrane. When the complex is bound to the test line, a burgundy colored band is visible, which confirms a reactive test result. Absence of a colored band in the test region indicates a non-reactive test result. To serve as a process control, a colored line in the control line area always changes from blue to red, indicating that the correct sample volume has been added and the membrane has been moistened.

### REAGENT COMPOSITION

25 sealed pouches each containing a test cassette, a dropper and a desiccant, 1 buffer and 1 package insert.

### MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection containers
- Centrifuge (for plasma only)
- Heparinized capillary tubes and dispensing bulb (for fingerstick whole blood only)
- Lancets (for fingerstick whole blood only)
- Timer

### REAGENT PREPARATION

The test is ready to use.

### STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. Do not freeze. Do not use beyond the expiration date.

### WARNINGS AND PRECAUTIONS

- For professional In Vitro diagnostic use only. Do not use after expiration date.
- This package insert must be read in full before performing the test. Failure to follow the insert will result in inaccurate test results.
- Do not use it if the tube/pouch is damaged or broken.
- Test is for single use only. Do not re-use under any circumstances.
- Handle all specimens as if they contain infectious agents. During the tests, follow the precautions against microbiological hazards and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Humidity and temperature can adversely affect results.
- Do not perform the test in a room with strong air flow.

### SPECIMEN COLLECTION AND STORAGE

- COVID-19 IgG/IgM Rapid Test Cassette can be performed using either whole blood, serum or plasma.
- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non hemolyzed specimens.
- Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimen may be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C. Whole blood collected by venepuncture should be stored at 2-8°C and the test must be carried out within 2 days. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be

tested immediately.

- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior for testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in accordance with local regulations for the transport of etiologic agents.

### TEST PROCEDURE

**Allow test cassette, specimen, buffer and/or controls to equilibrate to room temperature (15-30°C) prior testing.**

- Remove the test cassette from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
- Place the test device on a clean and level surface.


#### For serum or plasma specimens:

Using the 5 µL mini plastic dropper provided, draw a serum/plasma sample until the Specimen Line is crossed as shown in the following image, and then transfer the drawn serum/plasma sample into the sample well (S). Then add 2 drops (about 80 µL) of sample buffer to the buffer well (B) immediately. Avoid air bubbles.

**Note: Practice a few times prior to testing if you are unfamiliar with the mini dropper. For better precision, transfer the specimen by using a pipette capable of delivering 5 µL of volume.**

#### For whole blood specimen:

Hold the 5 µL mini plastic dropper vertically and transfer 1 drop of whole blood (about 10 µL) to the specimen well (S) of the test device. Then add 2 drops (about 80 µL) of sample buffer to the buffer well (B) immediately. Avoid air bubbles. Wait for colored line(s) to appear. The result should be read in 10 minutes. Positive results can be seen after 2 minutes. Do not interpret the result after 15 minutes.

### INTERPRETATION OF RESULTS

**Negative:** The colored line in the control line region (C) changes from blue to red. No line appears in the test line regions T1 or T2. The result is negative.

**IgM positive:** The colored line in control line region (C) changes from blue to red and a colored line appears in test line region T1. The result is anti-COVID-19 IgM positive.

**IgG positive:** The colored line in the control line region (C) changes from blue to red and a colored line appears in test line region T2. The result is anti-COVID-19 IgG positive.

**IgG and IgM positive:** The colored line in the control line region (C) changes from blue to red and two colored lines appear in test line regions T1 and T2. The result is anti-COVID-19 IgM and IgG positive.

**Invalid:** Control line is still completely or partially blue and fails to completely change from blue to red. Insufficient specimen volume or procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

### QUALITY CONTROL AND CALIBRATION

A procedural control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms a sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit. However, it is recommended to test positive and negative controls as a good laboratory practice to confirm the test procedure and to verify proper test performance.

### PERFORMANCE CHARACTERISTICS

#### Clinical Performance

The COVID-19 IgG/IgM Rapid Test has been evaluated with 113 blood samples obtained from patients exhibiting pneumonia or respiratory symptoms. The results were compared to RT-PCR or clinical diagnosis (including chest Computer Tomography and clinical signs etc.) of "Diagnosis and treatment of Novel Coronavirus pneumonia".

Regarding the IgM test, the result comparison to RT-PCR.

Method	RT-PCR		Total
	Positive	Negative	
COVID-19 IgG/IgM Rapid Test	87	0	87
Positive	12	14	26
Total	99	14	113

Regarding the IgG test, we counted the positive rate of the 36 out of 113 patients during the convalescence period.

Method	Number of patients during the convalescence period		Total
	Positive	Negative	
COVID-19 IgG/IgM Rapid Test	35	1	36
Positive	35	1	36
Total	35	1	36

The sensitivity of IgM test is 87.9% and specificity is 100% (14/14) compared to RT-PCR. The sensitivity of IgG test is 97.2% (35/36) during the convalescence period and specificity is 100% (14/14).

#### LIMITATIONS


- Use fresh samples whenever possible. Frozen and thawed samples (especially repeatedly) contain particles that can block the membrane. This slows the flow of reagents and can lead to high background color, making the interpretation of results difficult.
- Optimal assay performance requires strict adherence to the assay procedure described in this insert. Deviations can lead to incorrect results.
- A negative result for an individual indicates the absence of detectable anti-COVID-19 antibodies. However, a negative test result can occur if the quantity of the anti-COVID-19 antibodies present in the specimen is below the detection limits of the assay, or if the antibodies to be detected are not present during the disease stage in which a sample is collected.
- Some specimens containing unusually high titer of heterophile antibodies or rheumatoid factor may affect expected results.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.

#### WASTE MANAGEMENT

Please refer to local legal requirements.

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