DIAQUICK Influenza Ag Dipstick

An immunochromatographic lateral flow device for the qualitative detection of influenza type A and B in human nasal or throat swab samples

REF

Z09440CE

Content

- 25 tests individually packed (25x REF Z09440B)
- 25x 0.3 mL buffer tubes in a sealed zip lock pouch
- 25 sterilized swabs, individually packed
- 3 control swabs (A+, B+, negative)
- 1 working rack
- 1 package insert

For professional in vitro diagnostic use only.

GENERAL INFORMATION

| | |
|-------------------|---|
| Method | immunochromatographic antigen detection |
| Shelf life | 24 months from date of production |
| Storage | 2-30 °C |
| Sample | human nasal and throat swab samples |
| Results | after 10 minutes at room temperature |

DIAGNOSTIC EXPLANATION

Influenza illness classically presents with sudden onset of fever, chills, headache, myalgias, and a non-productive cough. Epidemics of influenza typically occur during winter months with estimated 114,000 hospitalizations¹ and 36,000 deaths² per year in the U.S. Influenza viruses can also cause pandemics, during which rates of illness and death from influenza-related complications can increase dramatically.

Patients who present with suspected influenza may benefit from treatment with an antiviral agent especially if given within the first 48 h after onset of the illness. It is important to rapidly distinguish influenza A from influenza B in order to allow physicians a choice to select the antiviral therapy. Moreover, it is important to determine if influenza A or B is causing the symptomatic disease in a particular institution (e.g. nursing home) or community, so that appropriate preventive actions can be taken for susceptible individuals. It is therefore important to not only rapidly determine whether influenza is present, but also which type of influenza virus is present.

TEST PRINCIPLE

The DIAQUICK Influenza Ag Dipstick uses monoclonal antibodies specific to Influenza type A and type B antigen for accurate determination of Influenza infection. Influenza antigens, type A and type B, in nasal or throat swab sample are allowed to react with the anti-Influenza type A and anti-Influenza type B monoclonal antibody-coupled gold conjugate followed by reaction with anti-Influenza type A and anti-Influenza type B monoclonal antibody in the test lines. When the nasal or throat specimen sample is infected with Influenza type A and/or type B, a visible line appears in the test region on the membrane. The DIAQUICK Influenza Ag Dipstick can also discriminate between Influenza type A and type B antigen.

REAGENT COMPOSITION

Influenza Ag Dipstick

Strip with the following components:

- Nitro-cellulose membrane strip with coated monoclonal antibodies against Influenza A and B in the test line regions and goat anti-mouse IgG in the control line region
- Conjugate pad with anti- Influenza A and B antibodies conjugated to colloidal gold particles
- sample and absorbant pad

Influenza Ag Extraction Buffer

- Biological buffer
- Sodium azide < 0.1 %

ADDITIONALLY REQUIRED MATERIAL

- General laboratory equipment

REAGENT PREPARATION

Before opening the pouch, the test must reach room temperature (18-25 °C). Once the pouch is opened, the test must be used immediately.

REAGENT STABILITY AND STORAGE

Conditions: The test must remain in the sealed aluminium pouch until use.

Storage: at 2-30 °C; Do not freeze!

Stability: up to the expiration date, when stored in the sealed pouch

Do not use the test, if the pouch is damaged.

SAMPLE COLLECTION AND PREPARATION

- **Nasal swab specimens:** Take out the nasal swab provided. Tilt the head of the patient backwards to the angle of 70° and gently insert the nasal swab into the nostril which shows the most secretion. Very gently rotate and push the swab until resistance is met at level of the turbinate. Gently rotate the swab against the nasal wall a few times.
- **Throat swab specimens:** To collect a throat swab specimen, vigorously rub a rayon throat swab on both tonsillar surfaces and the posterior pharynx. Remove swab from mouth and insert, tip down, into paper wrapper.

SAMPLE STABILITY AND STORAGE


Specimen should be tested as soon as possible after collection.

Do not use any kind of transport media to store or transport sample.

Specimens may be stored at 2-8 °C, or at room temperature (18-25 °C), in a clean, dry, closed container for up to 8 h prior to testing.

TEST PROCEDURE

1. Bring reagents and samples to room temperature (18-25 °C). The sensitivity of the test may be reduced at low temperatures.
2. Open the sealed buffer pouch, take out one buffer tube, then reseal the pouch via zip lock.
3. Peel off the protective cap of the buffer tube and insert the swab sample to be tested. Swirl the swab at least five times while pressing the head against the inside of the test tube.
4. Press the swab head against the inside of the buffer tube as you remove it. Dispose of the used swab in accordance with your biohazard waste disposal protocol.
5. Remove the test from the sealed pouch.
6. Place the test strip into the test tube with the arrows on the test strip pointing down. Do not handle or move the test strip until the test is complete and ready for reading.
7. Read the results **after 10 minutes**. False positive results can appear if the test is read after >15 minutes.

READING AND INTERPRETATION

NEGATIVE RESULTS

A coloured band is visible only in the control region (C).

POSITIVE RESULTS

Influenza Type A positive: two coloured bands are visible in the test region 1 (T1, adjacent to sample pad) and control region (C).


Influenza Type B positive: two coloured bands are visible in the test region 2 (T2, centre of the membrane) and control region (C).

INVALID RESULTS

If there is no coloured line in the control region (C), the result is invalid. This is due to deterioration of the test device or improper test procedure. Repeat the test with a new test strip.

QUALITY CONTROL

A procedural control is included in the test. A visible coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and a correct procedural technique.

It is recommended to test the positive and negative control swabs, which are supplied with this kit, as a good laboratory practice to confirm the test procedure and to verify proper test performance.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- Do not use the test strip after expiry date.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. The test strip should be discarded in a proper biohazard container after testing.
- The test strip is sensitive to humidity as well as to heat. So, it's very important to only take the test device from the sealed pouch shortly before it is used.
- Do not open the test foil pouch until ready to perform the test.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

LIMITATIONS

- The DIAQUICK Influenza Ag Dipstick is designed as a primary screening test of Influenza virus type A and type B infection.
- This test strip can provide a fast and easy way to get a result, but does not completely exclude the possibility of false positive or false negative results caused by various factors. So, the test results must be evaluated in conjunction with other clinical data available to the physician.
- Specimens that are collected after this phase may contain antigen titres below the reagent's sensitivity threshold. So a negative test results does not exclude infection with Influenza type A or type B.

PERFORMANCE CHARACTERISTICS

| Viral Culture / RT-PCR | | DIAQUICK Influenza Ag | | Commercial Influenza test | | |
|------------------------|-----------|-----------------------|-----|---------------------------|-----|-----|
| | | POS | NEG | POS | NEG | |
| POS | Type A | 75 | 69 | 6 | 69 | 6 |
| | Type B | 25 | 23 | 2 | 23 | 2 |
| | Sub Total | 100 | 92 | 8 | 92 | 8 |
| NEG | | 100 | 1 | 99 | 1 | 99 |
| TOTAL | | 200 | 93 | 107 | 93 | 107 |
| Sensitivity | | 92% (92/100) | | 92% (92/100) | | |
| Specificity | | 99% (99/100) | | 99% (99/100) | | |
| Accuracy | | 95,5% (191/200) | | 95,5% (191/200) | | |

Analytical specificity and cross reactivity

The DIAQUICK Influenza Ag Dipstick was evaluated with a total of 28 bacterial and viral isolates.

The studies demonstrated that there is no cross reactivity with the 28 tested bacterial and viral panels on the DIAQUICK Influenza Ag Dipstick.

Analytical sensitivity

The analytical sensitivity was established using a total of 6 strains of influenza viruses: 3 influenza A and 3 influenza B.

Precision

- 1) Within run performance test was determined by one analyst with 3 different lots
- 2) Between run performance test was determined by three analysts with 3 different lots containing different concentrations of each influenza type. There was no variation between different analysts.

LITERATURE

1. Comparison of Lateral-flow Immunoassay and Enzyme Immunoassay with Viral Culture for Rapid Detection of Influenza Virus in Nasal Wash Specimens from Children. Journal of Microbiology, Vol. 41, No. 5, may 2003, p. 2132-2134
2. Comparison of Two Nested PCR, Cell Culture, and Antigen Detection for the Diagnosis of Upper Respiratory Tract Infections due to Influenza Viruses. Journal of Medical Virology 59: 215-220 (1999)
3. Simultaneous Detection and Typing of Influenza Viruses A and B by a Nested Reverse Transcription-PCR: Comparison to Virus Isolation and Antigen Detection by Immunofluorescence and Optical Immunoassay (FLU OIA). Journal of Microbiology, Vol. 39, No. 1, Jan 2001, P. 134-138

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