

DIAQUICK COVID-19 AG CASSETTE

GENERAL INFORMATION

- The 2019 new coronavirus pneumonia (NCP), or “COVID-19”, was discovered for the 2019 Wuhan Viral Pneumonia case in China and was named by the World Health Organization on January 12, 2020.
- Coronaviruses are enveloped RNA viruses that are distributed broadly among humans, other mammals and birds and cause respiratory, enteric, hepatic and neurological diseases.


DIAQUICK COVID-19 AG FEATURES

- In-vitro immunochromatographic assay
- Rapid, qualitative detection of nucleocapsid protein antigen from SARS-CoV-2
- Swab from the nasopharynx

- **Fast results after 15 minutes**
- **Simple, time-saving procedure**
- **High sensitivity and specificity**



TEST PROCEDURE AND INTERPRETATION



ORDER INFORMATION

REF	Name	Specimen Type	Tests/Kit
Z20401CE	DIAQUICK COVID-19 Ag Cassette	swab	20